

**Contact [defibrillators.com.au](http://defibrillators.com.au) (1800 637 827) immediately if any problems are found**

**AED Monthly Inspection Form – Keep with AED**

Scan for  
AED information  
and forms



**Facility:**\_\_\_\_\_ **AED Serial No. (optional):**\_\_\_\_\_

**AED Location:** \_\_\_\_\_

[illegible]

**\*\*Required equipment and forms (must be with AED at all times):**

- 2 sets of adult chest pads - 1 set installed in the AED tightly, 1 back-up set - pads must be replaced if package has been opened
- 1 set of infant/child chest pads (except for Lifepak CR2 AEDs, which do not need separate infant/child pads)
- AED Monthly Inspection Form
- Resuscitation kit (2 pairs latex-free disposable gloves, CPR mask, pediatric CPR mask in pouch, shears, razor)